

Automated External Defibrillator

DefiMax AED



**European quality,
worldwide recognition**

Over 20 years of experience, thousands of patient monitors and defibrillators sold in Europe, Asia, Australia, Africa and South America.

• Emtel worldwide

Automatic external defibrillator (AED) dedicated to be used in public access areas as shopping malls, airports, railway stations, workplaces, sports facilities, etc.

Compliant with directions of the European Resuscitation Council. Clear voice messages and graphic instructions lead a rescuer easily through the whole cardiopulmonary resuscitation (CPR).

Easiness of use, intuitive interface and biphasic technology ensure the highest reliability of the device.

Technical parameters DefiMax AED



Polska

CE 2274

GENERAL PARAMETERS

Power supply	12V battery
Safety class	BF
Safety requirements	MDD: 93/42/EEC, 2007/47/EC; EN 60601-1, EN 60601-1-2, EN 60601-2-4
Weight	3,1 kg without electrodes
Dimensions	275 x 330 x 120 mm

MODES

AED mode

DEFIBRILLATOR MODULE

Electrodes type	single patient use
Impulse type	biphasic with patient's impedance compensation
Energy range	Adults up to 300 J Children up to 70 J

ADDITIONAL FUNCTIONALITY

CPR metronome
Self-test
Visual battery check

ACCESSORIES

Single patient use electrodes for children
Single patient use electrodes for adults
Battery Pack
Transport bag



A. Mickiewicza 66, 41-807 Zabrze, POLAND
phone: +48 32 271 90 13 fax: +48 32 271 57 27
medical@emtel.com.pl www.emtel.com.pl



Project co-financed by the European Union from the European Regional Development Fund.

Emtel reserves the right to make changes to this sheet and the product at any time without notice.